



DIAGNOSTICS

The Spirometry Source

Astra[★]300™ Features and Benefits

- | | |
|--|--|
| ■ Multifunction testing capability | ■ Tailor your testing modality as needed |
| ■ Program flexibility | ■ Choice of custom or pre-set programs |
| ■ Bright LED touchscreen display | ■ Facilitates patient data entry. Easy-to-read data, comments and interpretation |
| ■ Flow-trigger, patient-initiated testing | ■ Enhanced patient cooperation |
| ■ AstraPro Software included | ■ Download data, perform testing and communicate with your EHR system |
| ■ Measures over 40 parameters of FVC, VC and MVV | ■ Choose only the values needed for diagnosis |
| ■ Optional Bluetooth to PC | ■ Speedy wireless interface for interactive testing |
| ■ Internal temperature sensor | ■ Automatic BTPS correction |
| ■ Choice of 14 Predicted Authors | ■ Race specific equations/Ethnic correction |
| ■ Direct USB PC and printer interface | ■ Fast and easy documentation |
| ■ Full size report printouts | ■ Necessary for reimbursement |
| ■ Firmware upgrades at no charge | ■ Maintain up-to-date standards and features |
| ■ Test quality assurance prompts | ■ Insures reliable test results |
| ■ Automatic interpretation | ■ Aids in rapid diagnosis |
| ■ AccuFlo Digital Turbine | ■ Accurate and precise measurements |
| ■ Automated maneuver acceptability | ■ Assures repeatability |
| ■ Attractively priced | ■ Superior investment to income ratio |
| ■ No proprietary disposables | ■ Low cost per test, lower cost of ownership |
| ■ Audible end-of-test alert | ■ Assures accurate results |